# WHICH TABLE IS RIGHT FOR YOU?

# Surgical C-Arm Table - 840

Three-axis motorized positioning plus X-Y free-float



## SPECIFICATIONS:

Dimensions: 104" | x 26" w (260 x 66 cm) with OR accessory rails Tabletop: 101" | x 24" w (252 x 61 cm) Catheter Tray Extension (removable): 18" | x 24" w (45.7 x 61 cm) non-radiolucent Radiolucent A rea:

51" I x 24" w (128 x 61 cm) unobstructed head to accessory rail 20" I x 17.5" w (50.8 x 44 cm) additional, between accessory rails Radiolucent A rea, Total Length: 71" I (178 cm) Tabletop Material: Carbon fiber with integral head section Mattress: Seamless, 2" thick (5 cm) Attenuation:

Tabletop: .7 mm Aluminum equivalence Mattress: .5 mm Aluminum equivalence

#### Motions:

Height Adjustable: 33.5" to 43.5" (84 to 109 cm) Tabletop X Motion (head-to-toe float): 20" (51 cm) Tabletop Y Motion (side-to-side float): 10" (25.4 cm) Trendelenburg: 0° to 20° Reverse Trendelenburg: 0° to 20°

Lateral Roll: 0° to 20°

Controls:

H and Control: height, lateral roll and Trendelenburg motions Foot Control: activates all motions (has release for float) Push Button: releases X-Y Tabletop motions

#### W heels:

H ead End – 3" (7.62 cm) swivel casters, integral locking system on base Foot End – 5" (12.7 cm) swivel casters with central locking A ccessory Rails: Standard OR accessory rails 49.5" x 1.12" x.375" (124 x 2.86 x.95 cm) mounted near foot end of table

Patient Restraints: Three body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg)

based on line voltage Power: 115 VAC or 230 VAC and Battery. Battery automatically charges

when table is plugged in. An additional battery can be charged with optional wall-mounted charger. Standard Accessory:

– Catheter Tray Extension (removable)

- Bellows Skirt (disposable, 250/roll) and dispenser

Warranty: Two year

Certifications: ETL and cETL Listed to UL-60601-1° , CAN/CSA° C22.2 N  $\alpha$ : 601-1-M90, EN 60601-1 and IEC 60601-1-2, CE MARKED

 BI058-840
 Table, Surgical C-Arm, Contour, 115V

 BI058-840-10
 Table, Surgical C-Arm, Rectangle, 115V

 BI058-845
 Table, Surgical C-Arm, Contour, 230V

 BI058-845-10
 Table, Surgical C-Arm, Rectangle, 230V

# Surgical C-Arm Table - 846

Affordably priced float-top table



## SPECIFICATIONS:

Dimensions: 104" l x 26" w (260 x 66 cm) with OR accessory rails Tabletop: 101" l x 24" w (252 x 61 cm) Radiolucent A rea:

S1" I x 24" w (128 x 61 cm) unobstructed head to accessory rail 20" I x 17.5" w (50.8 x 44 cm) additional, between accessory rails Radiolucent A rea: Total Length: 71" I (178 cm) Tabletop Material: Carbon Fiber with integral head section Mattress: Seamless, 2" thick (5 cm) Attenuation:

Tabletop: .7 mm Aluminum equivalence Mattress: .5 mm Aluminum equivalence

Motions: Height: Adjustable from 33.5" to 43.5" (84 to 109 cm) Tabletop X Motion (head-to-toe float): 20" (51 cm) Tabletop Y Motion (side-to-side float): 10" (25.4 cm)

Controls: Hand Control: activates height Push Button: releases X-Y Tabletop motions

W heels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base Foot End – 5" (12.7 cm) swivel casters with central locking A ccessory Rails: Standard OR accessory rails 49.5" x 1.12" x .375" (124 x 2.86 x .95 cm) mounted near foot end of table Patient Restraints: Three body straps Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating. Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage Power: 115 VAC or 230 VAC Warranty: Two year Certifications: ETL and cETL Listed to UL-60601-1", CAN/CSA° C22.2

No.: 601-1-M90 and EN 60601-1, CE MARKED

Pain Management C-Arm Table - 870

Ideal for pain care applications



## SPECIFICATIONS:

Dimensions: 93" | x 26" w (233 x 66 cm) with OR accessory rails Tabletop: 90" | x 24" w (228.6 x 61 cm) Radiolucent Area: 60" | x 24" w (152.4 x 61 cm) unobstructed head to accessory rail 4" | x 17.5" w (10.1 x 44 cm) additional, between accessory rails Radiolucent Area, Total Length: 64" | (160 cm) Tabletop Material: Carbon fiber with integral head section Mattress: Seamless, 2" thick (5 cm) Attenuation: Tabletop: .7 mm Aluminum equivalence Mattress: .5 mm Aluminum equivalence

Motions:

Height Adjustable: 30" to 40" (75 to 100 cm) Trendelenburg: 0° to 20° Reverse Trendelenburg: 0° to 20° Lateral Roll: 0° to 20°

Controls: Hand Control: activates height, lateral roll and Trendelenburg motions Foot Control (optional): activates height, lateral roll and Trendelenburg motions

W heels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base Foot End – 5" (12.7 cm) swivel casters with central locking Accessory Rails: Standard OR accessory rails 30" x 1.12" x .375" (75 x 2.86 x .95 cm) mounted near foot end of table Patient Restraints: Three body straps Finish: Powder coat Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating. Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage Power: 115 VAC or 230 VAC Warranty: **Two year** Certifications: ETL and cETL Listed to UL-60601-1", CAN/CSA° C22.2 Na: 601-1-M90, EN 60601-1 and IEC 60601-1-2, CE MARKED

 BI058-846
 Table, Surgical C-Arm, Contour, 115V

 BI058-846-10
 Table, Surgical C-Arm, Rectangle, 115V

 BI058-847
 Table, Surgical C-Arm, Contour, 230V

 BI058-847-10
 Table, Surgical C-Arm, Rectangle, 230V

 BI058-870
 Table, Surgical C-Arm, Contour, 115V

 BI058-870-10
 Table, Surgical C-Arm, Rectangle, 115V

 BI058-875
 Table, Surgical C-Arm, Contour, 230V

 BI058-875-10
 Table, Surgical C-Arm, Rectangle, 230V

# Urology C-Arm Table - 800

## Configure to your exact specifications



#### SPECIFICATIONS:

Dimensions: 84.5" I x 25.5" w (214.6 x 64.8 cm) with OR accessory rails Tabletop: 53" | x 21" w (134.6 x 53.3 cm); 84" | (213.4 cm) with extension Lightweight Extension (for patient transfer only, not for imaging): 26.5" | x 21" w (67.3 x 53.3 cm)

Radiolucent Extension: 31" I x 25.5" w (78.7 x 64.8 cm) with OR accessory rails

Radiolucent Area: 54" | x 17.5" w (137.2 x 44.5 cm) Tabletop Material: Three Fiberesin<sup>™</sup> sections – each with individual 2" mattress

Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table Mattress: Individual, seamless, 2" thick (5 cm), one for each section

#### Motions:

Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm) Tabletop X Motion (head-to-toe): 20" (51 cm)

- Tabletop Y Motion (side-to-side): 10" (25.4 cm)
- Trendelenburg: 0° to 20° Reverse Trendelenburg: 0° to 20°
- Lateral Roll: 0° to 20°

Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

Controls

Hand Control: activates all motions Foot Control: activates all motions

#### Optional

X - Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.

X - Ray Grid: Removable X - Ray grid, 103 line, 8:1 ratio W heels:

Foot End – 3" (7.62 cm) locking swivel casters

Head End - 5" (12.7 cm) swivel casters with central locking A ccessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension) Patient Restraints: Two body straps, two ankle straps Finish: Powder coat

Patient Capacity:

I n urologic position: 500 lb (227 kg) I n use with radiolucent extension: 400 lb (181.4)

Weight tested to four times the patient load rating. Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.

Warranty: Two year Certifications: ETL and cETL listed to UL-60601-1, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1 and IEC- 60601-1-2:2001-09, CE MARKED

BI058-800 Table, Urology C-Arm - 800, 115 VAC BI058-805 Table, Urology C-Arm - 800, 230 VAC

# Brachytherapy C-Arm Table - 810

For image-guided seed implantations



### SPECIFICATIONS:

Dimensions: 84.5" | x 25.5" w (214.6 x 64.8 cm) with OR accessory rails Tabletop: 53" | x 21" w (134.6 x 53.3 cm); 84" | (213.4 cm) with extension Lightweight Extension (For patient transfer only, not for imaging.): 26.5" | x 21" w (67.3 x 53.3 cm) Radiolucent Extension: 31" I x 25.5" w (78.7 x 64.8 cm) with OR

accessory rails Radiolucent Area: 54" | x17.5" w (137.2 x 44.5 cm)

Tabletop Material: Three Fiberesin<sup>™</sup> sections – each with individual 2" mattress

Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table Mattress: Individual, seamless, 2" thick (5 cm), one for each section

Motions:

Ions: Height A djustable: 34" to 43.5" (86.4 to 110.5 cm) Tabletop X Motion (head-to-toe): 20" (51 cm) Tabletop Y Motion (side-to-side): 10" (25.4 cm) Trendelenburg: 0° to 20° Reverse Trendelenburg: 0° to 20° Lateral Roll: 0° to 20°

Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

Controls:

Hand Control: activates all motions Foot Control: activates all motions

Optional:

X -Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table. X - Ray Grid: Removable X - Ray grid, 103 line, 8:1 ratio

W heels: Foot End – 3" (7.62 cm) locking swivel casters Head End – 5" (12.7 cm) swive casters with central locking A ccessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension) Patient Restraints: Two body straps, two ankle straps

#### Finish: Powder coat Patient Capacity:

In urologic position: 500 lb (227 kg)

In use with radiolucent extension: 400 lb (181.4) Weight tested to four times the patient load rating. Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with

optional wall-mounted charger. Warranty: Two

Certifications: ETL and cETL listed to UL-60601-1, CAN/CSA® C22.2 No: 601-1-M90, EN 60601-1 and IEC- 60601-1-2:2001-09, CE MARKED

BI058-810 Table, Brachytherapy - 810, 115 VAC BI058-815 Table, Brachytherapy - 810, 230 VAC 0

# 3D Imaging C-Arm Table - 820

Designed to function with a 3D C-Arm



#### SPECIFICATIONS:

Dimensions: 103.75" | x 26" w (264 x 66 cm) with OR accessory rails Tabletop: 103.75" | x 24" w (264 x 61 cm) includes head and leg extension Radiolucent Area:

3D I maging A rea: 22.5" | x 24" w (57.1 x 61 cm) Radiolucent Extension (Fiberesin™ , non-3D): 29" | x 22.5" w (73.6 x 57.1 cm)

Total Radiolucent Area (Combined): 51.5" | (131 cm) Head-End Extension (non viewable): 18.2" | x 24" w (46 x 51 cm) Base Clearance to Floor: 1.5" h (4 cm) Tabletop Material: Carbon fiber

Mattress: Vinyl covered, seamless, 2" thick (5 cm)

Motions: Height Adjustable: 34" h to 44" h (86 to 112 cm) Tabletop X Motion (head to toe): 10" (25 cm) Tabletop Y Motion (side to side): 10" (25 cm) Trendelenburg: 0° to 20°

Reverse Trendelenburg: 0° to 20° Lateral Roll: 0° to 20°

Controls:

Hand Control: activates all motions Foot Control: activates all motions

## W heels:

Head End: 3" (7.6 cm) swivel casters, integral locking system on base Foot End: 5" (12.7 cm) swivel casters with central locking A ccessory Rails: Standard OR accessory rails, 1.125" x .375" (2.86 x .95 cm). Head-End 30" I (75 cm); Foot-End 7.75" I (20 cm) Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating. Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger. Weight: 542 lb (246 kg) Shipping Weight: 684 lb (310 kg) Warranty: Two y

Certifications: ETL and cETL listed to UL-60601-1, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1 and IEC- 60601-1-2:2001-09, CE MARKED

BI058-820 Table, 3D Imaging C-Arm, 115 VAC

BI058-825 Table, 3D Imaging C-Arm, 230 VAC