Surgical Case Carts



Standard Features:

- Fully-welded construction
- Sound-proofed tops and doors for quiet operation
- 8.5" x 11" card holder
- Adjustable perforated stainless steel roll out shelves
- Exterior push handles for better control and more interior space
- Non marring wrap-around bumper
- Positive lock door latches
- Sealed doors for sterility control
- Stainless steel casters
- Large 6" diameter, easy rolling, cart-washable casters with sealed bearings

Options:

- Fully extended perforated roll out shelves
- Top extensions helves
- Directional locks
- Total locks
- Silicon door gasket for complete sterility control
- Door safety seal tab to ensure door close through transport



Model	Case Dimension			OA Dimension			Adjustable Shelves	Card Holder
	Width	Depth	Height	Width	Depth	Height		
FHCMMCC1	26	24.25	47.25	31.75	28.375	56.375	2	1
FHCMMCC2	28.5	28.5	47.25	34.25	32.625	56.375	2	1
FHCMMCC3	40.25	24.25	30.25	46	28.375	39.375	1	1
FHCMMCC4	40.25	24.25	47.25	46	28.375	56.375	2	1
FHCMMCC5	36	25	47.25	41.75	29.125	56.375	2	1
FHCMMCC6	23.625	25	30.5	29.375	29.125	39.625	1	1
FHCMMCC7	35.5	24	44.5	41.25	28.125	53.625	2	1
FHCMMCC8	23	24.25	28.5	28.75	28.375	37.625	1	1
FHCMMCC9	23.625	25	46.25	29.375	29.125	55.375	2	1
FHCMMCC12	36	25	31	41.75	29.125	40	1	1
FHCMMCC14	46.5	25.25	31.25	52.25	29.375	40.25	2	1







Surgical Case Carts

Open Case Carts:

- Overall size: 40.375"W x 28.625"D x 41.125" H
- Case dimensions: 36"W x 25"D x 30.125"H
- Open design
- Full wrap around non-marring bumper
- 1 stainless steel adjustable shelf w/ fixed stainless steel top and bottom
- Stainless steel guard rail
- Large 6" diameter, easy rolling, cart-washable casters with sealed bearings



- Wire grid on back side
- 8" casters
- Wire shelves
- Push handles



Custom sizes available upon request

Open Case Cart w/ Wire Shelves:

- Overall size: 60.375"W x 26.125"D x 67"H
- Heavily built distribution cart
- Perfect for storage, supplies, linen or distribution
- Full wrap around non-marring bumper
- 3 standard solid stainless steel adjustable shelves
- Large 6" diameter, easy rolling, cart-washable casters with sealed bearings

Options:

- Wire grid on back side
- 8" casters
- Wire shelves



Custom sizes available upon request