

IV Stands

Constructed with fully welded stainless steel and tip-resistant base, the FHC IV stands are built with durability and reliability in mind. The FHC1315 and FHC1415 series IV poles give customers the option between hand or foot operated adjustment, and the number of hooks to be included.

FHC1315 Series Standard Features:

- Fully welded stainless steel construction
- 2" swivel casters
- Tip resistant base
- 68" - 108" height adjustment
- Hand operated
- Five leg stand

FHC1315-2 - 2 hook design

FHC1315-4 - 4 hook design

FHC1315-6 - 6 hook design


FHC1415 Series Standard Features:

- Fully welded stainless steel construction
- 2" swivel casters
- Tip resistant base
- 68" - 108" height adjustment
- Foot operated
- Five leg stand

FHC1415-2 - 2 hook design

FHC1415-4 - 4 hook design

FHC1415-6 - 6 hook design



FHC1410 Series Standard Features:

- 2" swivel casters
- Tip resistant base
- 68" - 108" height adjustment
- Hand operated
- Four leg stand
- Knocked-down design for low-cost shipping

FHC1410-2 - 2 hook design

FHC1410-4 - 4 hook design

FHC1410-6 - 6 hook design



FHC1510 Series Standard Features:

- 2" swivel casters
- Tip resistant base
- 68" - 108" height adjustment
- Hand operated
- Five leg stand
- Knocked-down design for low-cost shipping

FHC1510-2 - 2 hook design

FHC1510-4 - 4 hook design

FHC1510-6 - 6 hook design



IV Stands

Stainless Steel Lift Assist IV Pole Standard Features:

- Eliminates back injuries
- Easily lifts heavy loads by simply turning the crank
- Height adjustable from 56" - 96" high
- 6 leg base offers supreme stability
- 3" scuff resistant casters provide unmatched maneuverability
- 2" square stainless steel tube base post
- 1 1/4" stainless steel extension post tube
- No electric motors or hydraulic pumps to be concerned about or replace
- Many accessories optional available (universal clamp, infusion pump frame, add-a-tray)

FHC251 - 4 hook design

FHC252 - 8 hook design

FHC253 - 28 hook design



Shown with optional:
FHC260 - add-a-tray
FHC259 - infusion pump frame
FHC257 - universal clamps

Stainless Steel 5-Leg Infusion Pump Stand Standard Features:

- 16" diameter space saving base decreases obstruction in the operating room
- Base weight of 21 pounds for unmatched stability
- 3" scuff resistant casters
- Stainless steel pole section & detachable top
- Height adjustable from 52" - 94"
- Soft touch friction knob
- Standard base color is charcoal
- Optional base colors available upon request: (Red, Green, Yellow, Blue)

FHC221 - 2 hook design

FHC222 - 4 hook design

FHC223 - 6 hook design

Infusion Pump Stand Color Options:

