

Medical Cart Accessories

Drawer Accessories:



Clear Plastic Dividers

two per kit. Includes dividers. Fits in all drawers except half-sized Unicart or Medrx drawers.

FDIV-TRAY



White Divider Box w/ Clear Lid

White plastic divider box with clear lid. Can be sealed with clear plastic seal.

21.75"W x 15.75"D x 2.5"H

FEMER-BOX



Cut & Clip Divider System

Score the plastic strips and snap them to create your own divider system. Fits in all drawers. Set includes 22" long white plastic strips & 8 clips. FDIV-CC2 includes 2" H strips/FDIV-CC4 includes 4" H strips

**FDIV-CC2/FDIV-CC4
FDIV-CLIP**



White Plastic Tray with Configured Bins

Can be easily reconfigured. Fits in all full-sized Unicart drawers.

22"W x 15"D x 2.5"H

FTT-1A



2.5" Deep Configurable Gray Plastic Divider Tray

2.5" deep, includes dividers. Fits 2 per drawer, and typically used in 3" and 4" deep drawers.

9.25"W x 15"D x 2.5"H

FWDC-2025



White Plastic Tray with Configured Bins

Can be easily reconfigured. Fits in all full-sized Unicart drawers.

22"W x 15"D x 2.5"H

FTT-1B



4.75" Deep Configurable Gray Plastic Divider Tray

4.75" deep, includes dividers. Fits 2 per drawer, and typically used in 6", 9", and 12" deep drawers.

9.25"W x 15"D x 2.5"H

FWDC-2050



Adjustable Divider System w/ Row of Ampule Holders

21.5"W x 15.25"D

FDIV-DRW12



Adjustable Divider System w/ Large Divider Sections

Fits in 6", 9" and 12" deep full-sized Unicart drawers. 21.5"W x 15.25"D

FDIV-DRW5



Adjustable Divider System w/ Large Divider Sections

21.5"W x 15.25"D

FDIV-DRW34



Ampule Holder Divider System

With white plastic trays. Fits in all full-sized Unicart drawers.

21.5"W x 15.25"D x 2.3"H

FDIVTRAY-AMP



Ampule Holder Divider System

With white plastic trays. Fits in all full-sized Unicart drawers.

21.5"W x 15.25"D x 2.3"H

FDIVTRAY-ANES

Medical Cart Accessories

Drawer Accessories: (Continued)



Standard Divider System
With white plastic tray. Fits in all full-sized Unicart drawers.
21.5"W x 15.25"D x 2.3"H

FDIV-TRAY-STD
also available in 5" high version
FDIV-TRAY-STD5



Front & Back Divider System for Unicart Drawers
NOT for Med Jr. Drawers. Includes 3 front to back panels to create up to 4 compartments in one drawer.

FDIV-UC6 - For 6"H Drawers
FDIV-UC9 - For 9"H Drawers



Complete Ampule Holder Divider System
With white plastic tray. Includes 3 rows of 17 ampule dividers.
21.5"W x 15.25"D x 2.3"H

FWMT-3




Clear Dividable Trays
Trays fit in Unicart drawers and Instant Access Series.

FER-TRAY3
FER-TRAY6
For dividing systems, purchase **FDIV-DRW12ER**, **FDIV-DRW3ER** or **FDIV-DRW5ER**



Slim Med-Jr. Ampule Holder Divider System
Full drawer tray with 12 ampule dividers, ampule tray, rail and 4 adjustable dividers.
Fits 1 per drawer.
16"W x 16"D x 2.3"H

FWMT-4




Sealable, Perforated Bags
Fits over divider trays. Package of 20.

FBAG-3
23"W x 17"D x 5.5"H
FBAG-6
23"W x 17"D x 7"H



Slim Med-Jr. Full Drawer Divider System
Full drawer tray with 2 rails and 4 adjustable dividers. Fits 1 per drawer.
16"W x 16"D x 2.3"H

FWMT-5



Long or Short, Clear Plastic Dividers
For use with **FDIV-DRW12**, **FDIV-DRW34**, **FDIVTRAY-STD**, **FDIVTRAY-AMP**, **FDIVTRAY-ANES**

FDIV-10L - Long dividers
FDIV-10S - Short dividers



Divider Kit for 3" Aluminum Unicart Drawers
Includes 1 front-to-back and one side-to-side divider.

FDIV-AU3



Divider Kit for 3" Aluminum Med-Jr. Drawers
Includes 1 front-to-back and one side-to-side divider.

FDIV-AJ3



Divider Kit for 6", 9", & 12" Aluminum Unicart Drawers
Includes 1 front-to-back and one side-to-side divider.

FDIV-AUG



Divider Kit for 6", 9", & 12" Aluminum Med-Jr. Drawers
Includes 1 front-to-back and one side-to-side divider.

FDIV-AJ6

Medical Cart Accessories

Cart Accessories:



Lightweight Backside Cardiac Board Mount

Lightweight plastic cardiac board and metal mounting brackets for back side of Steel Unicart. Please see FCBB-7 for Aluminum Unicarts, FCBB-PC345 for FPC343 & FPC345. 24"W x 1"D x 18"H

FCBB-6



Lightweight Side Cardiac Board Mount

Lightweight plastic cardiac board and side mounting brackets for back side of Aluminum Unicarts or Med Jr. 24"W x 1"D x 18"H

FCBB-8



Thin Solid Plastic Cardiac Board Backside Mount

Thin solid plastic cardiac board and brackets for Steel or Aluminum Unicarts. 21"W x 15.25"D x 0.5"H

FCBB-9



FER2000 Backside Cardiac Board Mount

Lightweight plastic cardiac board mounting bracket and board for back side of the FER2000 cart.

FCBB-ER



Non-Locking Sharps Container Mounting Bracket

Mounts to side of all Unicarts, Med Jr, Medtrix, and DMC. Designed for 5 qt. container.

FSB-2



Utility Hooks

Mounts to existing holes on either side or back of cart. 2 hooks per package.

FUH-2



Oxygen Tank Mount

Plastic oxygen tank brackets with adjustable Velcro straps. For all Unicarts or Med Jr carts. For FPC343 or FPC345.

Please use part # FOTB-PC345.

FOTB-4



Metal Oxygen Tank Bracket

Can be mounted to all Unicarts and Med Jr carts.

FOTB-5



Sharp Container Mounting Bracket

With adjustable velcro strap. Mounts to side of all Unicarts, Med Jr, FER-2000, Medtrix, & DMC.


FSB-1



Locking Sharps Container Mounting Bracket

Mounts to side of all Unicarts, Med Jr, Medtrix, & DMC. Designed for 5 quart container.

FSB-2L



Self-Adhesive Chrome Hooks

Set of two self-adhesive chrome hooks that easily mount on any cart.

FABH-2



IV Pole

Chrome IV pole and mounting bracket. Mounts on Unicart, Med Jr, Medtrix, and DMC. For FER2000 please use part # FIV-3

FIV-2

Medical Cart Accessories

Cart Accessories: (Continued)



Coat Rack with Hooks
Made from lightweight aluminum. Mounts on Unicart, Med Jr, Medtrx, or DMC carts.

FCR-1



Side Mount Bucket Holder
Side mounting bucket holder and towel bar. Holds (2) 165oz cast buckets. 17.5"W x 6.5"D x 8.8"H

FCAST-1



Glove Box Holder
The NEW glove box holders are available in 3 sizes to fit any facilities' needs.

FBH-1 - one box (10"H x 4"D x 6"W)
FBH-2 - two box (10"H x 4"D x 12"W)



Small Waste Basket
Small (6 qt.) plastic waste basket. Requires accessory rail. 9"W x 6.25"D x 8.8"H

FBSK-1



Horizontal Glove Box Holder
Horizontal glove box holder. Holds 3 standard sized glove boxes.

FBH-3H



Large Waste Basket
Large (20 qt.) plastic waste basket. Requires accessory rail. 13"W x 8.5"D x 15"H

FBSK-2



Glove & Mask Glove Box Holder
Made from lightweight aluminum. Adjusts to many sizes of boxes. Mounts to railing on cart.

FGMH-1



Tilt Bins
Clear plastic tilt bins with white exterior. Available in 3,4,5, or 6 bin units. Mounts to base rail systems (FBRS-2, FBRS-3, or FBRS-JR)

FTB-3
FTB-4
FTB-5
FTB-6



Utility Receptacles
Available in 3 sizes. Requires accessory rail.

FUR-1 - .85 liters
FUR-2 - 1.5 liters
FUR-3 - 2.8 liters



Chart Holder
Gray steel chart holder. Mounts to accessory rail. 14.25"W x 1.75"D x 10"H

FCHT-1



Cylindrical Utility Receptacle
4" opening, 11" tall

FUR-4L



Notebook Mount
Wire notebook holder. Mounts to accessory railing. 12.5"W x 2.5"D x 10"H

FWNH-1

Medical Cart Accessories

Cart Accessories: (Continued)



Extra-Large Wire Basket
Mounts to accessory railing.
15.5"W x 9.5"D x 5.5"H

FWSB-1



Large Wire Basket
Mounts to accessory railing.
15.5"W x 5.5"D x 3.5"H

FWB-1L



Small Wire Basket
Mounts to accessory railing.
8.5"W x 5.5"D x 3.5"H

FWB-1S



Large Clear Plastic Organizer
Includes clear slide-in dividers. Mounts to accessory railing.
18"W x 5"D x 3"H

FORG-1L



Flat Shelf
Mounts to the horizontal rail system. (FBRS-2). Aluminum carts purchase FSHLF-3 & FBRS-3.
24"W x 10"D x 3/4"H

FSHLF-2



Two-Tier Shelving Unit
Mounts above Steel Unicarts. Order FSHLF-110A for aluminum carts.
26.5"W x 15.5"D x 1"H


FSHLF-110 - Steel
FSHLF-110A - Aluminum



Small Clear Plastic Organizer
Includes clear slide-in dividers. Mounts to accessory railing.

9"W x 5"D x 3"H

FORG-1S



Tape Dispenser

White tape dispenser with suction cups to attach to the top of a cart.


FTD-1



Horizontal Rail System

For hanging accessories and mounting shelves on Unicarts. Includes 2 vertical and 3 horizontal rails.

FBRS-2 - Steel or Aluminum carts
FBRS-JR - Med Jr. carts



Wrap-Around Aluminum Railing

May be used to hang accessories.

FR-3 - Steel carts
FR-4 - Aluminum carts
FR-JR - Med Jr. carts



Swiveling Monitor Shelf

Full swivel monitor shelf capable of locking in place. Includes adjustable strap to secure items in place. Max weight capacity 25 lbs. 13"W x 17"D x 12"H

FER-SHLF - Steel cart
FER-SHLFA - Aluminum cart
FER-SHLF2000 - Instant Access cart



Steel Scope Holder

Our durable steel scope holders secure easily to the wall or our Unicarts. Available in 2 sizes. Locking Plexiglas doors also included for added security.

FSH-2 - 2 scope holder
FSH-4 - 4 scope holder

Medical Cart Accessories

Cart Accessories: (Continued)



Collapsible Side Shelf

For Steel Unicart & Med Jr. Hangs on either side of cart. Surface area: 17.5"W x 15"D

FSSG-2C



Suction Shelf

For Unicart or Med Jr. Fits majority brands of suction units. 18" W x 8"D x 3"

FSS-M

Suction Pump

AC/DC portable unit provides 550mm HG and 50 lpm flow at the pump

FCEV3



Catheter Holder

Catheter holder for Unicart or Med Jr. 7.5"W x 2"D x 23"H

FCH-3



Backside Mounted Shelving System

This shelving unit includes (2) vertical rails, (1) top horizontal rail and (2) suture shelves. Aluminum carts require purchase of FAMP-2 for installation of this item.

Each shelf: 21.5"W x 5.5"D x 3.25"H

FSUT-1



Accessory Back Rail for Steel Med-Jr

Accessory rail for back attachment of accessories to the Steel Med Jr only.

FCB-2



Outlet Power Strip

Hospital grade outlet strip. 15 ft. long cord with 6 outlets.

FOS-1



Accessory Mounting Panel

For Lightweight Aluminum Unicarts. Mounts to SIDE to hold additional accessories.

FAMP-1



Accessory Mounting Panel

For Lightweight Aluminum Unicarts. Mounts to BACK to hold additional accessories.

FAMP-2



Accessory Side Rail

For Steel Unicarts or Steel Med-Jr carts. Rails for mounting various accessories to side of cart.

FSR-1



Accessory Back Rail

For Steel Unicarts or Steel Med-Jr carts. Rails for mounting various accessories to back of cart.

FCB-1



Cord Wrap

Usually used with outlet strip sold separately.

FCW-1



Fluid/Equipment Storage Tray System

For Steel, Med Jr or Medtrx/Halftr carts. Includes 2 trays. Each tray measures 17.52"W x 4.04"D x 4.5"H

FT-1 - Steel cart

FTA-1 - Aluminum cart

Medical Cart Accessories


Cart Accessories: (Continued)



Plastic Cart Cover

For Unicarts. Includes 5 bags.
Dimensions: 33" x 23" x 42"

FBAG-COVER5



Self-Adhesive Label Holders

1"H x 3"W, package of 100. Holders are loaded from the side.

FLH-1



Flexible Magnetic Light

For use on steel carts. Requires 2 AA batteries. Uses highly efficient LED bulb (included).

FLT-1



Locking Security Box

For use of transportation of narcotics or permanently mounted inside drawer. Includes 2 keys.
15.6"W x 11.4"D x 2.5"H

FSD-3



Deep Drawer Security Box

With security lid to secure medication in any full size Steel or Aluminum Unicart drawer 9"H. Includes lock w/ 2 keys.

FNARC-BOX



Proximity Lock Management Software

Management software for electronic and proximity locking systems.

FLOCKMANAGER

Proximity Card (pictured)

FPROX-1



Push-Button Locking Security Box

For use in transportation of narcotics. Permanently mounted. Includes override lock with 2 keys.

FSD-3S



Full Drawer Locking Security Box

Key locking, full drawer security box. Permanently mounted into drawer. Includes 2 keys. 21.5"W x 15.5"D x 2.5"H

FSD-3F




Aluminum Cart Gate Lock

For use on Lightweight aluminum carts. Can be used for securing narcotics and liquids or as a visible inventory management system. Available in Key-locking or Breakaway locking and in 3" increments up to 30" long.

FGLB -*Breakaway gatelock

FGLK -*Key gatelock




Steel Cart Gate Lock

For use on Steel carts. Available in Key-locking or Breakaway locking and in 3" increments up to 30" long.

FGLB -*Breakaway gatelock

FGLK -*Key gatelock



Self-Adhesive Proximity Tag

Small, adhesive proximity tag. 1" diameter

FPROX-2



Plastic Breakaway Seals


Can be used with lever locking systems, breakaway gate locks and FEMER-BOX divider trays.

FPS-100

Medical Cart Accessories

Example Cart Configurations:

FTRL-43336P-RED



FTGEA-333369-DKB



FTGK-433369-DKB-ANES-PKG



FTGLA-9PEDS-PKGA