



FUTURE HEALTH CONCEPTS, INC.

specializing in new and remanufactured medical equipment

SS Series Scrub Sink Standard Features:

- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve, one per station
- Check stop strainers for ease of maintenance
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sink are equipped with a 24 volt plug-in transformer
- Splash screen dividers with double and triple station sinks
- Main support and Z-brackets (in-wall support optional)
- 304 stainless steel polished and ground to a #4 finish

ES Series Scrub Sink Standard Features:

- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve, one per station mounted under sink
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sinks are equipped with a 24 volt plug-in transformer
- Z brackets (see specs on rear panel for in wall-support)
- 304 stainless steel polished and ground to a #4 finish

SETTING THE INDUSTRY STANDARD

1211 E. 30th Street • Sanford, FL 32773 • 800.282.8644 • 407.322.3871 fax • www.fhcusa.com