

# Ultrasonic Cleaners

## Ultrasonic Cleaning Systems



System Specifications:		CSU1	CSU3
Tank Dimensions	W	9 3/8 (23.8)	15 1/2 (39.4)
	H	4 (10.2)	6 (15.2)
	D	5 3/8 (13.7)	9 (22.9)
External Dimensions	W	12 5/8 (32)	19 (48.3)
	H	8 3/4 (22.2)	13 1/2 (34.3)
	D	9 7/8 (25)	16 (40.6)
Tank Capacity	Gallons (Liters)	0.8 (3.4)	3 (12.7)

Tuttnauer's Clean & Simple Ultrasonic Cleaning Systems are manufactured with only the highest quality components. These machines provide 60 kHz of quiet cleaning power which is less aggressive to instruments while still removing all surface contamination. The contemporary design and quiet operation compliments any office environment. Clean & Simple Ultrasonic Systems are available in 2 sizes with a full selection of accessories to meet your specific needs. The 1 gallon system holds a standard selection of instruments or 2 half size cassettes. The 3 gallon system holds 3 full size cassettes. Both sizes are available with an optional heater.

### FEATURES:

- High efficiency 60 kHz of cavitation for a thorough and quieter cleaning.
- Durable vinyl clad steel housing... easy to clean and resists scratching.
- Stainless steel tank with cover.
- Full line of accessories available.

## Ultrasonic Cleaning Tablets

Clean & Simple is an ultrasonic / enzymatic solution dispensed in a unique tablet form. Drop 2 tablets into your ultrasonic tank already filled with tap water and the tablets will dissolve in about a minute. This is a non-ionic solution which is non-corrosive, non-toxic, biodegradable and non-caustic. Available in boxes of 64 tablets, Item No. CS0064 and 144 tablets, Item No. CS0144.

### FEATURES:

- Easy to use tablet form.
- Eliminates handling of heavy 1 gallon containers.
- No bulky storage.
- No pre-mixing required.
- Enzymatic solution effectively breaks down protein.



TUCS0144

# Ultrasonic Cleaners

The FHC Cavitator Ultrasonic Cleaners come in 5 different sizes, designed to fit your cleaning needs. These cleaners generate ultrasonic waves that create millions of tiny bubbles through the process of cavitation. The bubbles gently clean and polish any surface they come in contact with, penetrating every nook and cranny to dislodge trapped dirt and debris. Medical instruments, glassware, even jewelry are quickly restored to their original lustrous finish without expending any effort. The gentle cleaning action is ideal for delicate items, too.




## Specifications:

Model	2.1	4.6	11	5.5S	18
Capacity	2.1 quarts (2 liters)	4.6 quarts (4.3 liters)	11 quarts (10.4 liters)	5.5 gallons (20.8 liters)	18 gallons (68 liters)
Outside Dimensions L x W x H					
inches	10 3/8 x 6 3/8 x 5 11/16	10 3/8 x 6 3/8 x 9 1/8	12 3/4 x 10 1/2 x 9	21 1/2 x 13 1/2 x 11 1/4	23 3/4 x 17 3/4 x 31
cm	26.4 x 16 x 14.4	26.4 x 16.2 x 23.2	32.4 x 26.7 x 22.9	54.6 x 34.3 x 28.6	60.3 x 45.1 x 78.7
Tank Size L x W x H					
inches	9 3/8 x 5 3/8 x 2 5/16	9 3/8 x 5 3/8 x 6	11 3/4 x 9 1/2 x 6	19 1/2 x 11 1/2 x 6	21 3/4 x 15 3/4 x 12
cm	23.8 x 13.7 x 5.9	23.8 x 13.7 x 15.2	29.8 x 24.1 x 15.2	49.5 x 29.2 x 15.2	55.2 x 40 x 30.5
Power Requirements	105/120 VAC 50/60 Hz	105/120 VAC 50/60 Hz	105/120 VAC 50/60 Hz	105/120 or 220 VAC 50/60 Hz	105/120 or 220 VAC 50/60 Hz
Power	85 watts	85 watts	200 watts	200 watts	800 watts
Nominal Main Frequency	67 KHz	67 KHz	67 KHz	67 KHz	24 KHz
Standard Accessories			Drain hose	Drain hose, Cleaning basket	Drain hose, Cleaning basket Fill hose, #1159



## Covers

- ME1012 Stainless steel cleaning baskets (fits ME11)
- ME1025 Stainless steel cleaning baskets (fits either ME2.1 or ME4.6)
- ME1038 Stainless steel cleaning basket (fits ME5.5S)



## Cleaning Baskets

- ME1011 Stainless steel cleaning baskets (fits ME11)
- ME1055 Stainless steel cleaning baskets (fits ME5.5S)
- ME1058 Stainless steel cleaning basket (fits either ME2.1 or ME4.6)
- ME1059 Stainless steel cleaning basket (fits ME18)