



2 PIECES OF HOOK AND LOOP TAPE ON BOTTOM OF PAD

NOTES:

- 1 CORE IS A HIGH DENSITY, FIRE RETARDANT POLYURETHANE FOAM
- 2 THE COVER IS FLUID PROOF, ANTI-MICROBIAL, FIRE-RETARDANT, AND STAIN RESISTANT.
- 3 COLOR: TOP AND SIDES, ROYAL BLUE, BOTTOM, BLACK

ANY PHYSICAL, DIMENSIONAL OR PERFORMANCE SPECIFICATIONS SHALL NOT BE REVISED WITHOUT WRITTEN AUTHORIZATION FROM FHC ENGINEERING

THIRD ANGLE PROJECTION INCHES

CRITICAL TO QUALITY C

TITLE
PAD, STRETCHER, 27 X 75 X 4

TOLERANCES UNLESS OTHERWISE SPECIFIED
 ANGLES: ±1° .XX = ±.030
 FRACTIONS: ± 1/64 .XXX = ±.015
 .XXXX = ±.005

DRAWN BY REX	DATE 7/9/2014
CHANGED BY	ECO #
CHECKED BY	DATE
SCALE 1:1	SHEET 1 OF 1

DRAWING NUMBER
FHC7200 MATTRESS

CONFIDENTIAL & PROPRIETARY TO FUTURE HEALTH CONCEPTS NOT TO BE COPIED OR SHOWN TO ANYONE WITHOUT THE EXPRESS PERMISSION OF FUTURE HEALTH CONCEPTS

REV.	DESCRIPTION	CHANGED BY:	DATE
A	RELEASED FOR PRODUCTION	REX	8/9/2006
REVISIONS			