



# FUTURE HEALTH CONCEPTS

“Your medical equipment specialist since 1975.”

## Warmer Storage Information

Model #	Total Cubic Feet	# of Blankets*	# of Solution 1 Liter Bottle	# of ½ liter bottles
FHCSWC1518	1.40	1-2	9	-
FHCSWC1824	3.00	6-9	20	-
FHCSWC24	5.27	9-14	36	-
FHCSWC2460	9.73	30-40	66	-
FHCSWC36	8.63	24-30	72	-
FHCSWC60	13.05	50-65	89	-
FHCSWC72	21.00	65-85	143	-
FHCDWC24 Lower Chamber	11.93	35-50	-	-
FHCDWC2424 Lower Chamber	11.93	35-50	-	-
FHCDWC2460 Lower Chamber	4.80	14-22	-	-
FHCDWC36 Lower Chamber	8.96	24-30	-	-
FHCDWC60 Lower Chamber	6.27	12-16	-	-
FHCTWC24 Lower Chamber	6.4	12-16	-	-
30” Roll-out Basket	-	-	30	42
1824 Roll-out Basket	-	-	16	20

**Please contact Future Health Concepts for any questions**

\* Number of blankets will depend on thickness of blankets, positioning, and fold.

**Setting the Industry Standard**